




# DU@Glide

## The Successful Synergy to Get Max Performance in Dermatological Applications

The Revolutionary Dual Wavelength Technology (CO<sub>2</sub>+1540 nm)

DUOGlide stands out for: Aesthetic Medicine



**DEKA**  
Innate Ability

# DUOGLIDE

EVOLVING EXCELLENCE



SKIN  
PHOTOREJUVENATION

ACNE SCARS

FACIAL WRINKLES  
REDUCTION

PERIOULAR  
LIFTING

EYELIDE PRO

DEEP SCARS


DERMATOLOGICAL  
SURGERY

# Key Principles

**DUOglide**, combines **two wavelengths (CO<sub>2</sub> 10600 nm and 1540 nm)** to maximize the efficacy in dermatological applications. Thanks to over 30 years of experience, DEKA has introduced a new system and series of more ergonomic and highly performing accessories.

The new **DUOglide** system exerts a single action on the tissue, effectively **stimulating collagen and minimizing downtime**. It is also perfectly suited to treatment of the most delicate areas like the neck, décolletage, and the area around the eyes.

## Why Choosing DUOglide:



# DU@Glide

THE PERFECT MATCH

## 1540 nm: The Ideal Partner for a New Era of Laser Therapy with CO<sub>2</sub> Laser Systems

Continuous research to maximize the efficacy, comfort, and safety of the treatment for the patient has led to the selection of two wavelengths: One ablative (CO<sub>2</sub>) and one thermal (1540 nm).



The second wavelength of 1540 nm available in the new miniaturized scanning systems was specifically developed to implement a synergy with the source CO<sub>2</sub> laser.

Thanks to the special sequential emission on the individual DOT, the effects of the two wavelengths are amplified synergistically, resulting in a reduced energy dosage but with the same efficacy, thus only a short downtime.

### Deeper Thermal Action

The synergy of the two CO<sub>2</sub>+1540 nm wavelengths also achieves heating, adjacent and non-coagulative of the entire scanning area, and reaches high dermal depth - not possible with only ablative sources.

The thermal effect reaches a depth level that maximizes tissue stimulation action and therefore obtains an even more effective treatment with reduced healing times.



THE PERFECT MATCH

## CO<sub>2</sub>+1540nm: Ideal Combination for Real Results

### Special Attention to Healing Times

Thanks to the deep and even stimulation of the special emission, the synergy of the two wavelengths boosts cellular turnover for faster healing, a real boon for patient recovery.

### Boosting of the Shrinkage Effect

The special sequential type emission achieved inside every individual DOT also synergically boosts the tissue shrinking effect to remodel and tone lax tissue.

### Emission Flexibility

Based on the type of application you can modulate the most suitable emission sequence: CO<sub>2</sub> + 1540 nm to boost the shrinkage effect (i.e. wrinkles and laxity), or 1540 nm + CO<sub>2</sub> for greater thermal effect for tissue stimulation.



THE SCIENCE BEHIND

# Maximum Efficacy, Precision and Total Control

## PSD Technology

The PSD (Pulse Shape Design) technology makes it possible to choose among various impulse modalities to be able to manage selectively the vaporization depth and the thermal effect: S-Pulse, D-Pulse, H-Pulse, U-Pulse, CW. By varying the impulse modality on the on the area concerned you obtain different ablation and stimulation effects to meet the various clinical needs.



SmartPulse (SP)




DEKAPulse (DP)



HighPulse (HP)

## SmartStack


For real and precise vaporization depth control you can select from 5 SmartStack levels. As the SmartStack level increases a progressive increase of the shrinkage effect is obtained, making the treatment more effective and safe, thus reducing patient recovery time.




## Smart Scanning Modality

### Scattered Modality

Scattered Modality evenly scans inside the entire area and fades along the edges to prevent overlapping of DOTs to ensure greater evenness between the scanning area and the skin texture.



### smart·TRACK



SmartTrack is the algorithm specially developed to optimize the scanning path to minimize local temperature increase.





THE SCIENCE BEHIND

## Set up Complete with Scanning Systems

DUOGlide has a range equipped with new scanning systems, optimized in shape, weight, easier connection and size to maximize handling.

### **μ-Scan DOT**

Scanning system designed to ensure utmost ergonomics in fractionated or traditional resurfacing treatments (modifiable parameters: Size, stretching and scanning area shape).



### **μ-SCAR 3**




### **μ-Scar 3**

Scanning system designed for remodelling deep scars. Minimizes the risk of post-treatment hyper and hypopigmentation thanks to the smaller spot size that generates greater depth of the action with reduced energy.

### **Derma Scan**


Scanning system paired with handpieces with 4" and 7" focus designed for fast and tridimensional ablations.




### **Eyelide PRO**

is the new laser blepharoplasty procedure which can be performed with both DuoGlide and Glide systems. By means of the "HP" pulse and its ergonomic "Slimcut" (2") handpiece we can obtain a precise incision of the eyelid in an operating field in which, thanks to the coagulating power of the CO<sub>2</sub> laser, bleeding is virtually absent. The treated tissue is stimulated through the tightening effect – which is why it is more elastic – benefiting from a better cicatrization process and lesser downtime, with a natural-effect result.

# Clinical Results




Courtesy of  
G. Scarcella, M.D.,  
Verona , Italy



Courtesy of  
G. Scarcella, M.D.,  
Verona , Italy



Courtesy of  
Prof. P. Bonan, M.D.,  
Florence , Italy



Courtesy of  
Prof. P. Bonan, M.D.,  
Florence , Italy





## The Floor to Practitioners

“ I have been using CO<sub>2</sub> laser systems since the mid 1990s when, while still a student, I was fascinated by this source and I have been constantly following its evolution ever since. The application versatility of these devices, with pulse modularity and controlled scanning, are today enhanced by a new action mechanism in which the 1540 nm wavelength increases the deep and homogenous thermal stimulation. The synergy of the two wavelengths, one with an ablative component and the other thermal, fosters cellular turnover with minimum downtime.”

**Prof. Paolo Bonan, MD- Dermatologist**  
*ESLD Key Officer, EADV Laser Task Force,  
Adjunct Professor of Laser at Plastic Surgery University of Siena  
In Charge of Laser Cutaneous Cosmetic & Plastic Surgery Unit  
Villa Donatello Clinic, Florence (Italy)*

“ I use the DUOglide Laser System (CO<sub>2</sub> & 1540 nm) and I find that, in addition to being extremely easy to handle, safe, and easy to use, for sure the most interesting and stimulating aspect is the real and effective synergy of the use of the two wavelengths practically simultaneously with the same device. This allows for greater effectiveness in terms of results both for treating scars of various origins and types and for photorejuvenation treatments of different body areas.”

**Dr. Giuseppe Scarcella, MD- Dermatologist**  
*Verona (Italy)*

THE PERFECT MATCH

## DUOglide Strength

- Maximum efficacy thanks to the **specially sequential emission** of the double wavelength.
- **Reduced downtime.**
- **Uncompromising versatility.**
- Exclusive advanced technology solutions **(PSD and SmartStack)** to maximize results
- Extreme flexibility of parameter settings
- Maximum **control and precision** of use.
- Wide use of applications thanks to **complete range of miniaturized scanning systems microscanners set up.**
- **Multi-decade experience in** CO<sub>2</sub> laser production



**Only for DEKA Users:**

A great support in medical practice

**DEKA** Club

**Intuitive GUI (15.6" rolling touch screen display)**

**PSD Technology**

**Miniaturized scanning systems**

**SmartTrack**

**Multimedia database**

**Teleassistance**



Scan and see  
**DUOglide** in action

# Technical Specifications

## DUOglide - Suggested configuration in Dermatology and Aesthetic Medicine

<b>CO<sub>2</sub> laser</b>	
<b>Laser Type</b>	CO <sub>2</sub> RF - PSD®
<b>Wavelength</b>	10,6 µm
<b>Laser emission mode</b>	TEM <sub>00</sub>
<b>Emission modes</b>	CW - SP - DP - HP - UP
<b>Average power at handpiece output</b>	0.1 - 60 W max
<b>IR laser</b>	
<b>Wavelength</b>	1540 nm
<b>Laser emission mode</b>	Circular multimode
<b>Power</b>	10 W
<b>General features</b>	
<b>Internal database</b>	More than 100 factory stored protocols, upgradeable with USB. Possibility for the user to store a number unlimited custom protocols
<b>Control panel</b>	Color LCD Touch Screen
<b>Accessories</b>	µ-Scan Dot, µ-Scar3, Dermascan, wide range of handpieces
<b>Dimensions and weight</b>	137 (A) x 42 (L) x 54 (P) cm - 70 kg

\*In this catalogue only the technical features of dermatology applications are listed.

DUOglide is also available in gynaecological configuration. Please refer to the related brochure for technical features in gynaecology.

**DANGER** - Visible and invisible laser radiation.  
Avoid eye or skin exposure to direct or scattered radiation.  
Class 4 laser product.

This brochure is not intended for the market of USA.

CE  
0123

DUOglide

**DEKA**  
Innate Ability

Follow us on  



  
[www.dekalaser.com](http://www.dekalaser.com)

Dealer stamp



**DEKA M.E.L.A. s.r.l.**  
Via Baldanzese,17 - 50041 Calenzano (FI) - Italy  
Tel. +39 055 8874942 - Fax +39 055 8832884

### DEKA Innate Ability

A spin-off of the El.En. Group, DEKA is a world-class leader in the design and manufacture of lasers and light sources for applications in the medical field. DEKA markets its devices in more than 80 countries throughout an extensive network of international distributors as well as direct offices in Italy, France, Japan and the USA. Excellence is the hallmark of DEKA's experience and recognition garnered in the sphere of R&D in over thirty years of activity. Quality, innovation and technological excellence place DEKA and its products in a unique and distinguished position in the global arena. DEKA manufactures laser devices in compliance with the specifications of Directive 93/42/EEC and its quality assurance system is in accordance with the ISO 9001 and ISO 13485 standards.